



















### DISPOSABLE ELECTRODES

#### SNAP ELECTRODES



##### FOR ADULTS

ref (part n.) description	picture	suggested use	ref (part n.) description	picture	suggested use
<b>F9060</b> 48x50mm FOAM solid gel		ECG recording long term monitoring holter stress test rest test	<b>F9049/100</b> 34x51mm FOAM solid gel		ECG recording long term monitoring holter stress test rest test
<b>F9079/100</b> 36x40mm FOAM solid gel		ECG recording holter stress test rest test	<b>F9069RM</b> 42x45mm FOAM solid gel radio-transparent		ECG recording rest test
<b>F9089/100</b> 36x45mm FOAM solid gel		ECG recording holter stress test rest test	<b>F9047</b> 32x44mm FOAM liquid gel		ECG recording holter stress test rest test
<b>F2080</b> Ø55 mm NON WOVEN solid gel		ECG recording long term monitoring holter stress test rest test	<b>F9047/4FL</b> 28x44mm FOAM liquid gel 4mm socket		ECG recording holter stress test rest test
<b>F9047/4F</b> 34x57 mm FOAM solid gel		ECG recording holter stress test rest test	<b>F9067</b> 48x50mm FOAM liquid gel		ECG recording holter stress test rest test
<b>F9048</b> 34x57mm FOAM solid gel		ECG recording holter stress test rest test			



##### PAEDIATRIC

ref (part n.) description	picture	suggested use	ref (part n.) description	picture	suggested use
<b>F9060P</b> Ø40mm FOAM solid gel		ECG recording long term monitoring holter stress test rest test	<b>PG10C</b> Ø26mm FOAM solid gel		ECG recording rest test
<b>F2080P</b> Ø40mm NON WOVEN solid gel		ECG recording long term monitoring holter stress test rest test	<b>F9069PRM</b> 32x36mm oval FOAM solid gel radio-transparent		ECG recording rest test
<b>F9079P/100</b> 32x36mm FOAM solid gel		ECG recording holter stress test rest test	<b>F9069NRM</b> Ø26mm FOAM solid gel radio-transparent		ECG recording rest test
<b>F9089P/100</b> 32x36mm FOAM solid gel		ECG recording holter stress test rest test			

### TAB ELECTRODES

ref	dimensions	backing	gel	radioltransparent	picture
<b>FOR ADULTS</b>					
F3001ECG	23 x 34 mm	plastic	solid	YES	
<b>PAEDIATRIC</b>					
F3001NECG	13 x 34 mm	plastic	solid	YES	

### PREWIRED ELECTRODES

ref	dimensions	backing	lead length	connector	gel	radioltransparent	picture		
<b>PER ADULTI - FOR ADULTS</b>									
F9053	Ø 32 mm	FOAM	35 cm	4 mm socket	solid	NO			
F9053/10			10 cm			NO			
F9053RM			35 cm			YES			
F9053RM/80			80 cm			YES			
<b>PAEDIATRIC</b>									
F9053N	23 x 30 mm	FOAM	35 cm	4 mm socket	solid	NO			
F9053N/10	23 x 30 mm		10 cm			NO			
F9053PR	20 x 15 mm		35 cm			NO			
F9053PR/10	20 x 15 mm		10 cm			NO			
F9053NRM	23 x 30 mm		35 cm			YES			
F9053NRM/80	23 x 30 mm		80 cm			YES			
F9053PRRM	20 x 15 mm		10 cm			YES			
F9058N/F2	23 x 30 mm		FOAM			50 cm		2 mm socket	NO
F9058P	Ø 32 mm					50 cm		NO	
F9058P/60	Ø 32 mm					60 cm		NO	
F9058N	23 x 30 mm		NOV WOVEN			50 cm		1,5 mm DIN socket	NO
F9058N/60	23 x 30 mm					60 cm			NO
F9058PR	20 x 15 mm					50 cm			NO
F9058PRM	Ø 32 mm					50 cm			YES
F9058NRM	23 x 30 mm					50 cm			YES
F9058/90NRM	23 x 30 mm					90 cm			YES
F9058PRRM	20 x 15 mm	50 cm		YES					
F9057P	Ø 32 mm	50 cm		NO					
F9057PRM	Ø 32 mm	50 cm		YES					
F9057N	23 x 30 mm	50 cm		NO					
F9057NRM	23 x 30 mm	50 cm	YES						
F9057PR	20 x 15 mm	50 cm	NO						
F9057PRRM	20 x 15 mm	50 cm	YES						
F9057PQ/60	32 x 32 mm	60 cm	NO						
F9057NQ/60	22 x 22 mm	60 cm	NO						

## REUSABLE ELECTRODES

### PERIPHERAL CLAMPS AND PLATE ELECTRODES


ref (part n.)	description	suggested use	picture
<b>FOR ADULTS</b>			
<b>F9024SSC</b>	clamp electrode Silver/Silver Chloride	ECG recording	
<b>F9024NI</b>	clamp electrode Nickel plated	ECG recording	
<b>F9024OSSC</b>	clamp electrode Silver/Silver Chloride, for Obese	ECG recording	
<b>F9024ONI</b>	clamp electrode Nickel plated, for obese	ECG recording	
<b>PAEDIATRIC</b>			
<b>F9023SSC</b>	clamp electrode Silver/Silver Chloride	ECG recording	
<b>F9023NI</b>	clamp electrode Nickel plated	ECG recording	
<b>PERIPHERAL PLATE</b>			
<b>F9010SSC</b>	Silver/Silver Chloride for adults	ECG recording	
<b>F9010PSSC</b>	Silver/Silver Chloride for paediatric	ECG recording	
<b>F9011PG</b>	rubber belt for peripheral plate electrode- h=24mm, L=48cm, for adults		
<b>F9011P</b>	rubber belt for peripheral plate electrode- h=17mm, L=26cm, for paediatric		

### SUCTION CHEST ELECTRODES

ref (part n.)	description	suggested use	picture
<b>with rubber ring</b>			
<b>F9002SSC</b>	Ø 24 mm, Silver/Silver Chloride	ECG recording	
<b>F9002NI</b>	Ø 24 mm, Nickel plated	ECG recording	
<b>F9003SSC</b>	Ø 30 mm, Silver/Silver Chloride	ECG recording	
<b>F9003NI</b>	Ø 30 mm, Nickel plated	ECG recording	
<b>with metallic ring</b>			
<b>F9008SSC</b>	Ø 15 mm, Silver/Silver Chloride	ECG recording	
<b>F9008NI</b>	Ø 15 mm, Nickel plated	ECG recording	
<b>F9015SSC</b>	Ø 20 mm, Silver/Silver Chloride	ECG recording	
<b>F9009SSC</b>	Ø 24 mm, Silver/Silver Chloride	ECG recording	
<b>F9009NI</b>	Ø 24 mm, Nickel plated	ECG recording	
<b>F9016SSC</b>	Ø 30 mm, Silver/Silver Chloride	ECG recording	
<b>F9016NI</b>	Ø 30 mm, Nickel plated	ECG recording	
<b>spare parts</b>			
<b>F9005SSC</b>	spare metallic chest, Ø 15mm, Silver/Silver Chloride	ECG recording	
<b>F9005/20SSC</b>	spare metallic chest, Ø 20mm, Silver/Silver Chloride	ECG recording	
<b>F9005/24SSC</b>	spare metallic chest, Ø 24mm, Silver/Silver Chloride	ECG recording	
<b>F9006SSC</b>	spare metallic chest, Ø 30mm, Silver/Silver Chloride	ECG recording	
<b>F9007</b>	spare rubber ball chest for series: F9003, F9009, F9016	ECG recording	
<b>F9007P</b>	spare rubber ball chest for series: F9002, F9008, F9015	ECG recording	

## PATIENT CABLES





In order to meet every customers' demand, FIAB produces a very wide range of ECG patient cables. Patient cables can be provided with 4 mm plugs, snaps or pinches at one side; at the opposite side of the cable different connectors are available allowing the connection with the most part of devices present on the market.



PATIENT DERIVATIONS		
small openable manifold, 4mm plugs 	large openable manifold, snap 	small openable manifold, clamps 





### ECG DEVICE MANUFACTURER

 <p>ASPEL BIONET BTL CARDIOLINE CARDIOTEK CARDIORAPID DR.LEE</p>	 <p>BURDICK / QUINTON</p>	 <p>ESAOTE</p>
 <p>BOSCH CAMBRIDGE KENZ NIHON KOHDEN PICKER</p>	 <p>CARDIOLINE</p>	 <p>FUKUDA DENSHI</p>
 <p>BOSCH CUSTO MED GMBH HELLIGE SIEMENS</p>	 <p>CARDIOLINE KENZ</p>	 <p>FUKUDA DENSHI</p>
 <p>BURDICK / QUINTON SIEMENS</p>	 <p>CARDIOLINE CARDIETTE</p>	 <p>HELLIGE MARQUETTE</p>
 <p>BURDICK / QUINTON SIEMENS</p>	 <p>ESAOTE</p>	 <p>MORTARA</p>



WITH INSPECTABLE MANIFOLD			WITH MOULDED MANIFOLD
<p>F6000 series, with small manifold and 4mm plug connection</p> 	<p>F6000P series, with small manifold and pinch connection</p> 	<p>F5000c series, with large manifold and snap connection</p> 	<p>F8000 series, with small moulded manifold and 4mm plug connection</p> 

## ACCESSORIES AND CONNECTION CABLES


ref (part n.)	description	picture
<b>KITS OF 10 UNIPOLAR LEADS</b>		
<b>F5377/KIT</b>	4mm plug connectors, without protection	
<b>F5388/KIT</b>	4mm plug connectors, 4,7 Kohm protection	
<b>F5388R/KIT</b>	4mm plug connectors, 10 Kohm protection	
<b>F5376/KIT</b>	Pinch connectors, without protection	
<b>F5386/KIT</b>	Pinch connectors, 4,7 Kohm protection	
<b>F5386R/KIT</b>	Pinch connectors, 10 Kohm protection	
<b>F5399/KIT</b>	Snap connectors, without protection	

<b>ADAPTERS</b>		
<b>F9017W</b>	From 4mm socket to snap+alligator, white	
<b>F9017/4PW</b>	from 4mm socket to pich, white	
<b>F9019</b>	from 4mm socket to alligator (*)	
<b>PG922/4T</b>	from 4mm socket to snap (**)	
<b>PG922/4KIT</b>	from 4mm socket to snap , kit of 10 pieces (***)	

(\*) suffix colour key: R=red, N=black - (\*\*) suffix colour key: R=red, N=black, V=green, G=yellow, B=white - (\*\*\*) 10pcs: 6 white, 1 black, 1green, 1yellow

<b>SINGLE SNAP CONNECTION CABLES</b>		
<b>F9022S</b>	L=60cm, 2mm safety plug (*)	
<b>F9052S</b>	L=80cm, 2mm safety plug (*)	
<b>F9032S</b>	L=150cm, 2mm safety plug (*)	
<b>F9022D</b>	L=60cm, safety 1,5mm DIN socket (*)	

(\*) suffix colour key: R=red, N=black, A= Orange, V=green, G=yellow, W=white, M=brown, B=blue, l=purple

<b>KITS OF SNAP CONNECTION CABLES</b>		
<b>F9022D/KITX</b>	Kit of dnap connection cables, L=60cm, safety 1,5mm DIN socket	
<b>F9052D/KITX</b>	Kit of dnap connection cables, L=90cm, safety 1,5mm DIN socket	
<b>F9032D/KITX</b>	Kit of dnap connection cables, L=120cm, safety 1,5mm DIN socket	

F9022D/KITX and F9052D/KIT are available with 3,5,7 or 8 cables, please replace the "X" with the related number. For model F9032D/KITX 8 cables version is not available